

Safe

Fast

Easy

Cell Total RNA Isolation Kit
培养细胞 RNA 提取试剂盒



Foregene Co.,Ltd.

产品订购热线: 400-8566189

地址: 成都市温江区八一路北段18号三医创新中心A区D2栋

邮箱: info@foregene.com

网址: www.foregene.com



福际
生物

专注生命科学 · 让实验变得更简单

产品简介 INTRODUCTION

福际的培养细胞总RNA试剂盒采用本公司研制的离心柱和配方,可以高效率的从96、24、12、6孔板培养细胞中提取得到高纯度高质量的总RNA。

试剂盒提供高效的DNA-Cleaning Column,能轻松的让上清液和细胞裂解物分离并吸附除去基因组DNA,操作简便、省时;

该RNA-only Column能高效的结合RNA,搭配独特的配方,可以同时处理大量样品。



无需低温冰浴

无需氯仿

无需DNase I


订购信息 ORDERING INFORMATION

产品名称	规格	货号
培养细胞总RNA提取试剂盒 Cell Total RNA Isolation Kit	50T	RE-03111
	200T	RE-03113


下游应用 APPLICATION

培养细胞总RNA提取试剂盒提取得到的总RNA可用于各种下游分子实验,例如:cDNA合成,RT-PCR、Real Time PCR、Northern Blot、Dot Blot、体外翻译、芯片分析、PolyA筛选、分子克隆和RNase保护分析等


实验流程 EXPERIMENTAL PROCESS



实验图例 EXPERIMENTAL LEGEND



Cell Total RNA Isolation Kit处理上述不同数量的细胞,20μl体积洗脱,取2μl纯化的总RNA 1%琼脂糖凝胶电泳图谱



细胞类型: MKN-74
细胞数量: 10⁶、10⁵、10⁴

引用文献 REFERENCES

- ◆ Yuan, B., Peng, Q., Cheng, J. et al. Structure of the Ebola virus polymerase complex. *Nature* 610, 394–401 (2022). Doi: 10.1038/s41586-022-05271-2. **IF: 69.504 (Foregene #RE-03111)**
- ◆ Meng, Y., Qiu, L., Zeng, X. et al. Targeting CRL4 suppresses chemoresistant ovarian cancer growth by inducing mitophagy. *Sig Transduct Target Ther* 7, 388 (2022). Doi: 10.1038/s41392-022-01253-y. **IF: 38.104 (Foregene #RE-03111)**
- ◆ Ren, Y., Wang, A., Wu, D. et al. Dual inhibition of innate immunity and apoptosis by human cytomegalovirus protein UL37x1 enables efficient virus replication. *Nat Microbiol* 7, 1041–1053 (2022). Doi: 10.1038/s41564-022-01136-6. **IF: 30.964 (Foregene #RE-03111)**

相关产品 RELATED PRODUCTS

第一链cDNA合成预混体系Master Premix (For qPCR)

去基因组第一链cDNA合成试剂盒 (For qPCR) (With gDNase)

荧光定量PCR试剂盒-SYBR Green I 染料法

鼠尾直接PCR试剂盒—防PCR产物污染体系

用过都说好!

Make the experiment easier